



Reactor-Ready[™] Lab Reactors





Innovative, patented, low cost, reaction work stations for glass vessels from 100 ml to 5 litres

Reactor-Ready Lab Reactor

Replace multiple reaction set-ups with a single universal system that accepts vessels from 100 ml to 5 litres















Unique vessel clamp

- self centering
- quick release
- change vessels in minutes

Hose couplings

- quick release
- · connect in seconds
- chemically resistant PEEK

Individual manifolds easy connection easy drain down reduces stress on vessel sidearms



Simple stirrer alignment slide and glide

- set and lock

Patented stirrer coupling

- drop in
- no tools
- no fuss

Standard vessels

- 250 ml to 5 litres
- vacuum jacketed option

Process vessels

- mimic plant scale geometry
 - 100 ml to 5 litres
- vacuum jacketed option

Reactor-Ready Duo accepts two vessels and allows reactions in parallel or series






Change to Reactor-Ready

It improves your chemistry by saving space, time and money

Optional AVA **Control Software** allows you to log and control stirrers, circulators, balances, pumps, temperature sensors and other devices





Accepts all leading brands of overhead stirrer





Compatible with all leading brands of circulator or thermoregulator

Radleys Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom

Glass and

PTFE lids





support collar



Reactor-Ready

Swap reaction vessels in minutes, NOT hours

Reactor-Ready is designed as a universal reactor work station that can be used for different vessel sizes and different experiments. Reactor-Ready can easily and quickly be configured to suit the chemistry and scale you need for each project.

The beauty of Reactor-Ready is that this one reactor work station can replace many, saving you money and fumehood space.

Features include

- · Rapid vessel exchange with quick-release vessel clamp and hose couplings.
- Range of single and vacuum jacketed vessels from 100 ml to 5 litres.
- Range of process vessels to mimic larger scale plant or manufacturing reactors.
- · Accepts all leading brands of overhead stirrer and allows easy, tool-free adjustment.
- · Triple support stand features heavy duty stainless steel support rods for strength and stability.
- Self aligning stirrer coupling engages without the need for tools.
- Temperature range: -70 °C to +230 °C (jacket temperature)
- · Innovative hose manifolds allow easy thermofluid drain down.
- Wide range of accessory glassware including condensers, dropping funnels etc.
- Optional software allows you to log and control stirrers, circulators, balances, pumps, temperature sensors and other devices.

Patented vessel clamp allows rapid vessel exchange without tools

- Unique clamp design allows vessel to be supported and clamped by a single operator in less than a minute.
- · Use the same clamp to hold all vessel sizes; plus the vessels remain supported even when the clamp is unlocked.











Vessel support collar aids vessel loading

- Reactor-Ready uses an innovative support collar fitted to the neck of each vessel, which allows the vessels to be swapped in and out of the Reactor-Ready core like a cartridge.
- The stainless steel collar locks securely into the support plate on the framework and allows the vessel to remain supported, even when not clamped.
- Vessel exchange is tool free and can easily be undertaken by one person.

Unique self aligning stirrer coupling saves time when changing vessels

- Patented drop-in stirrer coupling engages easily without the need for tools.
- Simply lower the stirrer coupling to engage with the stirrer guide below.
- The stainless steel stirrer coupling will fit in any standard stirrer chuck (for 10mm diameter)
- Powerful stirring up to 500 rpm, or 800 rpm for short periods (subject to motor selected).





Choice of PTFE stirring shafts

- Reactor-Ready is offered with a choice of PTFE stirring shafts: anchor, retreat curve and turbine/propeller styles.
- PTFE anchor stirring shaft is included with all standard reaction vessel kits
- PTFE turbine stirring shaft is included with all process reaction vessel kits.

I-beam support allows the stirring motor to be positioned with ease

- I-beam allows the stirrer motor to be raised or lowered without needing to remount the stirrer.

 Adjustable perfets step college legated below the Library allow the upper to lower the stirrer.

 Adjustable perfets step college legated below the Library allow the upper to lower the stirrer.

 Adjustable perfets step college legated below the Library allow the upper to lower the stirrer.
- Adjustable safety stop collars located below the I-beam allow the user to lower the stirrer motor to a set position for ease of operation and safety.
- Reactor-Ready accepts all leading brands of motor or air-powered overhead stirrer with support rods up to 16 mm in diameter.





PTFE Stirrer Guide

- The Reactor-Ready's PTFE stirrer guide screws directly on to the glass lid and features a unique extended support which minimises wobble even at high speeds.
- The stirrer guide is chemically inert, gas-tight and suitable for work under vacuum.

Quick release PEEK hose couplings

- Operates over the temperature range -70 °C to +230 °C (jacket temperature).
- · Quick-release, screwthread design allows rapid exchange of vessels without the use of tools.
- Simply tighten by hand for a leak-tight seal.
- PEEK with Viton O-rings ensures all wetted parts are chemically resistant.
- Lightweight, metal-free design prevents chipping and minimises strain on the glass fitting.
- · All hoses have a 12 mm ID, maximising thermofluid flow and heating/cooling performance
- Low thermal conductivity of PEEK reduces frosting when cooling.

Hose manifolds allow easy thermofluid drain down

- The inlet and outlet hoses from the Reactor-Ready vessel connect to manifolds attached to the support rod. These manifolds serve two important functions:
- Support the hoses running from the vessel, thereby eliminating stress on the glass sidearms.
 Allows thermofluid drain-down, which provides a quick and clean method of removing the
- Allows thermofluid drain-down, which provides a quick and clean method of removing the thermofluid from the vessel jacket when changing vessels.
- By opening the venting valve on the (upper) outlet manifold, thermofluid drains from the vessel back to the thermoregulator under gravity.
- Opening the drain tap on the (lower) inlet manifold allows drainage to an external reservoir.





curve stirrer

er

propeller stirrer











Case study: Accelerating process optimisation in a busy chemical development lab

Researchers at a leading US pharmaceutical company invested in Reactor-Ready for use in their chemical development lab. The process R&D group, specialising in small molecule research, are using Reactor-Ready for optimising new or existing chemistry, identifying and eliminating impurities, and screening process variables for transfer to kilo-labs, pilot plant and manufacturing sites.


Commenting on their existing methods the Lead Scientist said, "We use round bottom flasks, but temperature control via mantles or oil baths is poor. Alternatively, when we use traditional jacketed reactor vessels, we need a different set up for every scale change. Even changing a reactor of the same size takes ages."

"So I liked Reactor-Ready the minute I saw how easy it was to change reactors. Reactor-Ready does it all. I get one stand which fits all sizes up to 5 litres. I get overhead stirring on all scales. I get vacuum jacketed vessels, which provide excellent temperature control with minimum heat loss, which is a real advantage at low temperatures, and I get the ability to rapidly change out my reactors."

"The first time I changed vessels, it took about 5 minutes but now I can change in less than 2 minutes. That includes removing and attaching the heating fluid hoses. Occasionally I lose a drop or two of heat transfer fluid but most times I lose none. The thing I like about the system is the control and speed. I get better control of my reactions and I can vary my work or scale immediately".

Radleys Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom





Reactor-Ready core system includes

- Stand with base, support rods, bosses and all clamps
- Glass reaction vessel lid
- · Vessel FEP O-ring seal and PTFE collar
- HDPE drip tray
- · Vessel-to-manifold insulated hoses with quick-release connections
- Lower inlet manifold + drain
- · Upper outlet manifold + vent
- PTFE stirrer guide
- Stainless steel drop-in stirrer guide coupling
- Two stainless steel vessel support collars
- · Locking rod for support rod bosses



Reaction vessel kits make ordering vessels simple and more cost effective

These cost effective vessel kits include the vessel and accessories you need to get started. The Pt100 probe and PTFE stirrer included in each kit are matched to the vessel size. Fach vessel kit includes:

- · Reaction vessel with PTFE stopcock
- Sidearm couplings
- Pt100 PTFE temperature probe with LEMO connector
- Probe adapter





stand is secure and stable

Reactor-Ready Duo

vacuum jacketed option

process and standard

geometry

Thermofluid manifolds

 choice of two manifold kits • single circulator manifolds double circulator manifolds

Supports two glass reaction vessels of different or same volume from 100 ml to 5 litres

Reactor-Ready Duo shares the same unique features as Reactor-Ready, except that it supports two separate jacketed glass reaction vessels. The system can be configured to operate with a single thermoregulator controlling the temperature of both vessels simultaneously, or with two thermoregulators controlling the temperature of each



Patented stirrer coupling drop in

set and lock

for each vessel

- no tools
- no fuss

Hose couplings

- quick release
- connect in seconds
- chemically resistant PEEK

Strong framework

- heavy duty base
- solid s/steel support rods
 - · large drip tray

Why Reactor-Ready Duo?

All the benefits of Reactor-Ready with two vessels in parallel or series

Features include

- · Rapid exchange of both vessels independently, with quick-release vessel clamp and hose couplings.
- Choice of manifold kits allow two vessels to run from a single thermoregulator or two separate thermoregulators.
- System accepts two overhead stirrers which can be operated and moved independently.

Applications

- Parallel synthesis or reaction optimisation: use the same or different size vessels and vary stirring speed, stirrer shape and temperature between vessels as appropriate.
- Two stage reaction: transfer reactant from one vessel to the other using vacuum or a pump.
- Single reaction vessel: using the second vessel as either a receiving or feed vessel (where reagents can be pre-heated or pre-cooled prior to addition).
- · Use optional AVA Software to control fluid transfer between vessels.









Similar vessels with







How to order Reactor-Ready Duo...







4. Select the overhead stirrers.



5. If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware. then select from the accessory list

6. If you need automation, add the AVA Software and data hub.



Reactor-Ready Duo core system

Easy-to-order Reactor-Ready Duo core system includes all the framework components required for operation

Features include

- · Robust heavy-duty framework features solid stainless steel support rods.
- Two separate quick release clamps support each of the vessels independently.
- · Large drip tray will retain vessel contents.
- Manifold and hose kits must be ordered separately.

Core system (RR121300) includes

- · Stand with base, support rods, bosses and all clamps
- · 2 x glass reaction vessel lids
- 2 x vessel FEP O-ring seals and PTFE centering collars
- 2 x PTFE stirrer guides and couplings
- · HDPE drip tray
- · 2 x stainless steel stirrer drop-in couplings
- · 2 x stainless steel vessel support collars
- · locking rod for support rod bosses



Reactor-Ready Duo manifolds and hoses

A choice of manifold kits for use with one or two thermoregulators

Double thermoregulator manifold and hose kit (RR121310)

This double thermoregulator kit is designed to work with two separate thermoregulators controlling the temperature of each vessel independently. With this set-up both vessels can operate at different temperatures. This is the recommended kit

- R/H upper return manifold with vent valve
- L/H upper return manifold with vent valve
- R/H lower supply manifold with drain valve
- L/H lower supply manifold with drain valve
- 4 x quick-release insulated hoses, 60 cm
- 4 x sealing caps for hoses



Single circulator manifold and hose kit (RR121305)

This single circulator kit is designed to work with one thermoregulator controlling the temperature of both vessels simultaneously. With this set-up both vessels would operate at the same temperature. The manifolds are supplied with blanking plugs should you wish to run only one vessel.

- · central upper return manifold with vent valve
- · central lower supply manifold with drain valve
- 4 x quick-release insulated hoses, 70 cm
- 4 x sealing caps for hoses
- 2 x stainless steel blanking plugs





Hei-Torque overhead stirrers

The powerful Hei-Torque range can accomplish the most demanding mixing tasks whilst providing the highest safety

Stirrer models include

- Powerful stirring from 10 to 2000 rpm (Max. recommended speed with Reactor-Ready is 500 rpm).
- Available with basic (Expert) or advanced (Ultimate) features, and a choice of up to 40, 200 or 400 Ncm torque
- Intuitive control panel for easy operation, safe start and stop to avoid accidental start-up.
- Ultimate range features: set speed limit, adjustable acceleration, set a torque limit.
- Unique Quick-chuck for immediate and convenient 'one-hand' impeller changes, without tools,

	Hei-Torque Core	Hei-Torque Expert 200	Hei-Torque Expert 400	Hei-Torque Ultimate 200	Hei-Torque Ultimate 400
Power rating, motor input/output (W)	105/75	120/80	150/90	120/80	150/90
Number of gears	1	1	2	1	2
Speed range (rpm)	20 to 2000	10 to 2000	10 to 400 and 200 to 2000	10 to 2000	10 to 400 and 200 to 2000
Speed indicator	digital monochrome 2"	digital monochrome 2.4"	digital monochrome 2.4"	digital colour 3.2"	digital colour 3.2"
Speed control	electronic	electronic	electronic	electronic	electronic
Torque, maximum (Ncm)	40	200	400	200	400
Torque indicator (Ncm)	symbol	symbol	symbol	precise value	precise value
Overheat protection	automatic cut-out	automatic cut-out	automatic cut-out	automatic cut-out	automatic cut-out
Viscosity up to (mPa s)	10,000	100,000	250,000	100,000	250,000
Stirring cap. H ₂ O max. (L)	25	50	100	50	100
Shaft diameter, max. (mm)	10.5	10.5	10.5	10.5	10.5
Dimensions (W x H x D) (mm)	70 x 281.5 x 195	86 x 350 x 247	93 x 350 x 247	86 x 350 x 247	93 x 350 x 247
Weight (kg)	2.3	5.1	5.3	5.1	5.3
Set speed limit, Torque limit, Adjustable acceleration	-	-	-	Yes	Yes
Analogue/digital interface	-	-	-	USB and RS232	USB and RS232


For Reactor-Ready the maximum speed is 500 rpm for continuous use for 800 rpm for short periods

Thermoregulators/circulators

Which thermoregulator should you choose?

- · Radleys are a Huber distribution and service partner, so we recommend the Huber range of thermoregulators, circulators, chillers and waterbaths. However, Reactor-Ready is compatible with all leading brands of thermoregulator including Julabo, Lauda, Polyscience, Haake and many others,
- Also available are a range of hoses, adapters, temperature probes and thermofluid for all applications.

Description	Temperature range °C	Heating power (kW)	Cooling power at 0 °C (kW)	Application Examples		
Huber Ministat 230 -45 °C to 200 °C	45 °C to 200 °C	2.0	0.38	250 ml	Cool to -10 °C in <1 hr, then heat to 130 °C in <30 mins	
	-45 C t0 200 C			1 litre	Control at 0 °C, then heat to 100 °C in 20 mins	
Huber Petite Fleur	-40 °C to 200 °C	1.5	0.45	500 ml	Cool to -20 °C in approx 30 mins, then heat to 150 °C in <30 mins	
	-40 C to 200 C			1 litre	Heat to 150 °C, then cool to 0 °C at > 2.5 K/min	
Huber	Huber Unistat Tango 45 °C to 250 °C	3.0	0.7	1 litre	Cool to -30 °C in <1 hr, then heat to 180 °C in <1 hr	
Unistat Tango				2 litre	Heat to 180 °C, then cool to 20 °C in approx 1 hr	
Huber Unistat 405	-45 °C to 250 °C	3.0	1.0	2 litre	Cool to -20 °C in <30 mins, then heat to 100 °C in <30 mins	
				5 litre	Control at 20 °C, then cool to -20 °C in approx 70 mins	





Huber Unistat 405



Huber Petite Fleur

Reaction vessels

Single jacketed, vacuum jacketed and process styles

Features include

- · Precision engineered borosilicate glass vessels can be used with Reactor-Ready and Reactor-Ready Duo.
- DN100 ground and polished flat flange combines with the FEP-encapsulated O-ring and our unique self centering PTFE collar to offer a leak-tight seal between the vessel and lid.
- Convenient and cost-effective vessel kits include stirrer, temperature probe and adapter.

Vessels feature zero dead space PTFE bottom outlet valve

- Exclusive stopcock design prevents over-tightening that can cause glass breakage.
- · Piston design allows quick and safe removal, eliminating the cleaning and assembly problems associated with other stopcock styles.
- Glass-filled PTFE piston minimises leaking or weeping across the full temperature range (-70 °C to +230 °C)
- Interchangeable between vessels with a full range of replacement parts and FEP O-rings.
- The large 15 mm bore minimises sample hold-up and is ideal for slurries (25 mm available on request).

Glass vessel lids

- DN100 5 neck glass lid fits all vessel sizes.
- DN100 lid with 5 x Rodaviss, grease-free, safety joints. B24 centre, B34, B29, B24 and B19 side joints.
- Rodaviss glass sockets will accept any glassware with traditional A, B or C length conical cones.

Standard Vessels















250 ml vacuum











500 ml single





iacketed vessel

500 ml vacuum

iacketed vessel





500 ml single jacketed

3000 ml single

PTFE vessel lids

DN100 5 neck PTFE lids fit all vessel sizes.

 Lids accept a range of removable Rodaviss PTFE sockets and compression fittings.

· More robust and less likely to break than glass.

• Ideal for use with analytical probes including Blaze.

 Rodaviss PTFE sockets will accept any glassware with traditional A, B or C length conical cones.



Process Vessels - mimic the geometry of plant scale reactors

Feature a 1.25 to 1 ratio of internal height to internal diameter, with dish-shaped bottom.







Note: 3 litre process vessels are not shown



100 ml vacuum jacketed







250 ml single jacketed














Reaction vessel kits

Reaction vessel kits make ordering vessels and accessories simple and much more cost effective






Glassware and accessories

Wide choice of standard glassware to suit your application

Glassware features Rodaviss joints and is compatible with Reactor-Ready lids









Reducing

adapter R29



adapter B29

adapter B29



Reflux divider head

B29 + GI 18



Reflux divider head

B29 + GI 18 + temp-

probe port B14



inlet adapter, B24







Glass cone

stopper, B29



Other styles and sizes of glassware are available, please see full listing.

Radleys Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom

How to order Reactor-Ready

- 1. Select the Reactor-Ready core: RR121000.
- 2. Choose the vessel kit or kits you require by selecting the vessel volume and style.
- Select the overhead stirrer you require.
- 4. If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list opposite.
- 5. If you need automation, add the AVA Software and Data Hub.

Reactor-Ready core system

- RR121000 Reactor-Ready core system includes:
 - Stand with base, drip tray, support rods, bosses and all clamps
 - · Vessel to framework insulated hoses with quick-release connections
 - Upper and lower drain-down manifolds
 - PTFE stirrer guide and drop-in coupling
 - Glass reaction vessel lid with Rodaviss sealing caps
 - Vessel FEP O-ring seal and PTFE centering collar
 - 2 x Stainless steel vessel support collars

Important: Hoses to connect the Reactor-Ready framework to your thermoregulator are NOT included. You will require hoses with an M24x1.5 female connection or appropriate adapters. See opposite for hoses, adapters and suitable thermoregulators.

Standard Reaction Vessel Kits

Kits include: standard vessel, PTFE anchor stirrer paddle, PTFE Pt100 temp probe (with LEMO plug for connection to the thermoregulator) and PTFE temp probe adapter.

RR120010 250ml Single Jacketed Vessel Kit RR120015 250ml Vacuum Jacketed Vessel Kit RR120020 500ml Single Jacketed Vessel Kit RR120025 500ml Vacuum Jacketed Vessel Kit RR120030 1000ml Single Jacketed Vessel Kit RR120035 1000ml Vacuum Jacketed Vessel Kit RR120040 2000ml Single Jacketed Vessel Kit RR120045 2000ml Vacuum Jacketed Vessel Kit RR120042 3000ml Single Jacketed Vessel Kit RR120047 3000ml Vacuum Jacketed Vessel Kit



Kit includes Anchor Stirre

Process Reaction Vessel Kits

Kits include: process vessel, PTFE turbine stirrer paddle, PTFE Pt100 temp probe (with LEMO plug for connection to the thermoregulator) and PTFE temp probe adapter.

RR120105 100ml Single Jacketed Process Vessel Kit RR120107 100ml Vacuum Jacketed Process Vessel Kit 250ml Single Jacketed Process Vessel Kit RR120110 RR120115 250ml Vacuum Jacketed Process Vessel Kit 500ml Single Jacketed Process Vessel Kit RR120120 RR120125 500ml Vacuum Jacketed Process Vessel Kit RR120130 1000ml Single Jacketed Process Vessel Kit RR120135 1000ml Vacuum Jacketed Process Vessel Kit RR120140 2000ml Single Jacketed Process Vessel Kit RR120145 2000ml Vacuum Jacketed Process Vessel Kit RR120142 3000ml Single Jacketed Process Vessel Kit RR120147 3000ml Vacuum Jacketed Process Vessel Kit RR120150 5000ml Single Jacketed Process Vessel Kit RR120155 5000ml Vacuum Jacketed Process Vessel Kit



Kit includes

How to order Reactor-Ready Duo

- 1. Select the Reactor-Ready Duo core: RR121300.
- 2. Choose manifold and hose kit for one or two circulators.
- 3. Choose the vessel kit or kits you require by selecting the vessel volume and style.
- 4. Select the overhead stirrers you require.
- 5. If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list opposite.
- 6. If you need automation, add the AVA Control Software and Data Hub.

Reactor-Ready Duo core system

RR121300 Reactor-Ready Duo core system includes:

- Duo stand with base, drip tray, support rods, bosses and all clamps
- 2 x PTFE stirrer guides and drop-in couplings
- 2 x Glass reaction vessel lids with Rodaviss sealing caps
- 2 x Vessel FEP O-ring seal and PTFE centering collars
- 2 x Stainless steel vessel support collars

Important: Hoses to connect the Reactor-Ready Duo framework to your thermoregulator are NOT included. You will require hoses with an M24x1.5 female connection or appropriate adapters. See opposite for hoses, adapters and suitable thermoregulators.

Double Circulator Manifold and Hose Kit

RR121310

- Double Circulator Manifold and Hose Kit includes:
- R/H Upper Return Manifold with Vent Valve
- L/H Upper Return Manifold with Vent Valve
- · R/H Lower Supply Manifold with Drain Valve L/H Lower Supply Manifold with Drain Valve
- 4 x Quick-Release Insulated Hoses, 60cm
- 4 x Sealing Caps for Hoses

Single Circulator Manifold and Hose Kit

- Single Circulator Manifold and Hose Kit includes:
- Central Upper Return Manifold with Vent Valve
- Central Lower Supply Manifold with Drain Valve
- 4 x Quick-Release Insulated Hoses, 70cm
- 4 x Sealing Caps for Hoses
- 4 x Stainless Steel Blanking Plugs

Hei-Torque Overhead Stirrers

501-60410-00 Hei-Torque Core 501-62010-00 Hei-Torque 200 Expert

501-62020-00 Hei-Torque 200 Ultimate (USB/RS232)

501-64010-00 Hei Torque 400 Expert

501-64020-00 Hei-Torque 400 Ultimate (USB/RS232)

General Maintenance Accessories

Reactor-Ready Maintenance Kit (includes items listed below) Oil Wipe for Support Rods, Pack 10

Stirrer Guide Bearing Seal B24 10mn RR121063

RR121064 Stirrer Guide Ferrule and Shaft Seal B24 10mm RR121071 23mm Viton O-Rings for Quick Release Hose Coupling, Pack 2

RR121072 28mm Viton O-Rings for Quick Release Hose Coupling, Pack 2

RR121090 FEP Encapsulated Silicone Vessel O-Ring Replacement 15mm Piston Top O-Ring RR166100 Replacement 15mm Piston Middle O-Ring RR166102

Reactor-Ready Duo reaction vessel compatibility guide

The following combinations of reaction vessel sizes can be used in the Reactor-Ready Duo.

	100ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml	
Standard Vessel, Single Jacketed								
250 ml		•	•	•	•	•		
500 ml		•	•	•	•	•		
1000 ml		•	•	•	•	•		
2000 ml		•	•	•	•	•		
5000 ml		•	•	•	•	•		
	Process Vessel, Single Jacketed							
100 ml	•	•	•	•	•	•	•	
250 ml	•	•	•	•	•	•	•	
500 ml	•	•	•	•	•	•	•	
1000 ml	•	•	•	•	•	•	•	
2000 ml	•	•	•	•	•	•	•	
3000 ml	•	•	•	•	•	•	n/a	
5000 ml	•	•	•	•	•	n/a	n/a	

	100ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml
	Standard Vessel, Vacuum Jacketed						
250 ml		•	•	•	•	•	
500 ml		•	•	•	•	•	
1000 ml		•	•	•	•	•	
2000 ml		•	•	•	•	•	
5000 ml		•	•	•	•	n/a	
		Proce	ss Vessel, \	Vacuum Jac	keted		
100 ml	•	•	•	•	•	•	•
250 ml	•	•	•	•	•	•	•
500 ml	•	•	•	•	•	•	•
1000 ml	•	•	•	•	•	•	•
2000 ml	•	•	•	•	•	•	n/a
3000 ml	•	•	•	•	•	n/a	n/a
5000 ml	•	•	•	•	n/a	n/a	n/a

Hubor Thorm	oregulators (Circulators)					
HB2015.0005.0						
HB2016.0005.01 Ministat 240 Pilot ONE Controller						
HB1030.0001.0						
HB1000.0037.0 HB1002.0045.0	3					
	ator Hoses, Adapters and Thermofluid					
	with PTFE Inner					
HB9325	M24x1.5 Insulated Hose NW12 100cm -60°C to +260°C					
HB9326 HB9327	M24x1.5 Insulated Hose NW12 150cm -60°C to +260°C M24x1.5 Insulated Hose NW12 200cm -60°C to +260°C					
HB9328	M24x1.5 Insulated Hose NW12 300cm -60°C to +260°C					
	with Metal Inner					
HB9274 HB9275	M24x1.5 Insulated Hose NW12 100cm -100°C to +350°C M24x1.5 Insulated Hose NW12 150cm -100°C to +350°C					
HB9276	M24x1.5 Insulated Hose NW12 200cm -100°C to +350°C					
HB9277	M24x1.5 Insulated Hose NW12 300cm -100°C to +350°C					
Adapters HB6724	M24x1.5 Female to M16x1 Male Adapter					
HB6945	M24x1.5 Male to M16x1 Female Adapter					
HB6723 HB6195	M24x1.5 Female to M30x1.5 Male Adapter M16x1 90° Adapter					
HB9256	M24x1.5 90° Adapter					
HB6089	M16x1 Nut					
HB12634 HB6086	M24x1.5 Nut M16x1 Hose Connector NW8					
HB9471	M24x1.5 Hose Connector NW12					
HB6091	M16x1 Isolating Ball Valve -20°C+140°C (max. 6 bar at +140°C)					
HB9236 Thermofluid	M24x1.5 Isolating Ball Valve -10°C+180°C (max. 6 bar at +180°C)					
HB6161	SilOil M20.195/235.20, 5 Litre -20°C to +195°C/235°C					
HB6162	SilOil M20.195/235.20, 10 Litre -20°C to +195°C/235°C					
HB6163 HB6164	SilOil M40.165/220.10, 5 Litre -40°C to +165°C/220°C SilOil M40.165/220.10, 10 Litre -40°C to +165°C/220°C					
HB6165	SilOil M60.115/200.05, 5 Litre -60°C to +115°C/200°C					
HB6166	SilOil M60.115/200.05, 10 Litre -60°C to +115°C/200°C					
HB6479	DW-Therm M90.200.02, 10 Litre -90°C to +200°C (Unistats only)					
Glassware Ad RR139005	Condenser Liebig 250mm B29 + GL14 + fittings					
RR139007	Condenser Double Surface 200mm B29 + GL14 + fittings					
RR139009	Condenser Jacket Coil Rodaviss B29 + GL14 + fittings					
RR139011 RR139013	Right Angle Adapter B29 + GL14 + PTFE fittings 50ml Graduated Dropping Funnel Equalising Arm B24					
RR139015	100ml Graduated Dropping Funnel Equalising Arm B24					
RR139017	250ml Graduated Dropping Funnel Equalising Arm B24					
RR139019 RR139021	500ml Graduated Dropping Funnel Equalising Arm B24 1000ml Graduated Dropping Funnel Equalising Arm B24					
RR139023	Swan Neck Adapter Angled B24 to Vertical B24					
RR139025 RR139027	Gas Purge Adapter B24 + GL14 + fittings Glass Cone Stopper B19					
RR139027 RR139029	Glass Cone Stopper B24					
RR139031	Glass Cone Stopper B29					
RR139033 RR139035	Glass Cone Stopper B34 Enlarging Adapter B24 to B29					
RR139037	Enlarging Adapter B29 to B34					
RR139039	Reducing Adapter B24 to B19					
RR139041 RR139043	Reducing Adapter B29 to B19 Reducing Adapter B29 to B24					
RR139045	Reducing Adapter B34 to B19					
RR139047	Reducing Adapter B34 to B24					
RR139049 RR139051	Reducing Adapter B34 to B29 Oil Bubbler Type 1 B24 + GL14 + fittings					
RR139053	Solid Additions Funnel B34 (with extended tube)					
RR139057 RR139063	Valved Liquid Inlet Adapter B24 to GL14 PTFE Hose Adapter B19 for 12mm ID Tubing					
RR139067	Boss Head					
RR139069	3-Prong Clamp					
RR139079 RR139081	PTFE Flexible Bellows B24 PTFE Flexible Bellows B29					
RR139083	PTFE Flexible Bellows B34					
RR139150	Reflux Divider Head B29 + GL18 + PTFE Fitting					
RR139152 RR139190	Reflux Divider Head B29 + GL18 + Temp Probe Port B14 Reflux Divider Kit					
RR139192	Reflux Divider Kit with Vapour Temp Kit					
RR139154 RR139324	250ml Receiver Flask 2 x B24 Socket, + B24 Cone and Tap 500ml Receiver Flask 3 x B24 Socket, + B24 Cone and Tap					
RR139324 RR139326	1000ml Receiver Flask 3 x B24 Socket, + B24 Cone and Tap					
PTFF Pt100	Temperature Probes and Adapters					
RR140250	Probe 8mm with LEMO Plug - 100ml, 250ml + 500ml Vessels					
RR141000	Probe 9mm with LEMO Plug - 1000ml, 2000ml + 3000ml Vessels					
RR145000 RR140010	Probe 9mm with LEMO Plug - 5000ml Vessel PTFE Temperature Probe Adapter - 8mm					
RR140015	PTFE Temperature Probe Adapter - 9mm					
Reaction Ves	sel Support Collar					
RR121048	Vessel Support Collar					

Glass Vessel Lid RR136000 DN100 Glass 5 Neck Lid - 1 x B19, 2 x B24, 1 x B29, 1 xB34 Includes Rodaviss sealing caps PTFE Vessel Lids DN100 PTFE 5 Neck Lid & 2 Blanked Ports. 1 x B19, 2 x B24, 1 x B29, 1 x B34 Includes PTFE sockets/blanking plugs and Rodaviss sealing caps DN100 PTFE Lid only with 7 Threaded Ports. 2 x 1/4" NPT, 5 x 3/4" NPT PTFE sockets are not included - select from below PTFE Vessel Lids - for Blaze Metrics Probe DN100 PTFE 5 Neck Blaze Lid & 2 Blanked Ports. 1 x B19, 2 x B24, 1 x B29, 1 x B34. Includes PTFE 14mm and 19mm Compression Fittings, PTFE sockets/blanking plugs and Rodaviss sealing caps DN100 PTFE Blaze Lid only with 7 Threaded Ports. 2 x 1/4" NPT, 5 x 3/4" NPT PTFE sockets are not included - select from below PTFE Rodaviss Sockets, Compression Fittings & Blanking Plugs for PTFE Lids Rodaviss B19 Socket Adapter (¾ NPT Rodaviss B24 Socket Adapter (¾ NPT Rodaviss B29 Socket Adapter (3/4" NPT Rodaviss B34 Socket Adapter (3/4" NPT RR136303 Rodaviss B45 Socket Adapter (3mm Bore Compression Fitting 1 RR136314 RR136315 RR136316 RR136317 8mm Bore Compression Fitting 12mm Bore Compression Fitting 14mm Bore Compression Fitting RR136319 16mm Bore Compression Fitting 3/4" NPT 19mm Bore Compression Fitting 3/4" NPT RR136320 PTFE Blanking Plug ¼" NPT PTFE Blanking Plug ¾" NPT RR136328 Standard Reaction Vessels Only 250ml Single Jacketed Vessel 250ml Vacuum Jacketed Vessel PP130500 500ml Single Jacketed Vessel RR130550 500ml Vacuum Jacketed Vesse 1000ml Single Jacketed Vessel RR131050 1000ml Vacuum Jacketed Vessel 2000ml Single Jacketed Vessel RR132050 2000ml Vacuum Jacketed Vessel 3000ml Single Jacketed Vessel 3000ml Vacuum Jacketed Vessel RR133000 PTFE Stirrer Shafts (for Standard Vessels) Anchor PTFE Stirrer Shaft - 250ml vessels Retreat Curve PTFE Stirrer Shaft - 250ml vessels
Turbine PTFE Stirrer Shaft - 250ml vessels Anchor PTFE Stirrer Shaft - 500ml vessels Retreat Curve PTFE Stirrer Shaft - 500ml vessels Turbine PTFE Stirrer Shaft - 500ml vessels RR150505 RR151000 RR151005 Anchor PTFE Stirrer Shaft - 1000ml vessels Retreat Curve PTFE Stirrer Shaft - 1000ml vessels Turbine PTFE Stirrer Shaft - 1000ml vessels RR152000 RR152005 Anchor PTFE Stirrer Shaft - 2000ml vessels Retreat Curve PTFE Stirrer Shaft - 2000ml vessels RR152010 Turbine PTFE Stirrer Shaft - 2000ml vessels RR153000 Anchor PTFE Stirrer Shaft - 3000ml vessels Retreat Curve PTFE Stirrer Shaft - 3000ml vessels Turbine PTFE Stirrer Shaft - 3000ml vessels Process Reaction Vessels Only 100ml Single Jacketed Process Vessel 100ml Vacuum Jacketed Process Vessel RR130152 250ml Single Jacketed Process Vessel RR130257 250ml Vacuum Jacketed Process Vessel RR130502 500ml Single Jacketed Process Vessel RR130552 500ml Vacuum Jacketed Process Vessel RR131002 1000ml Single Jacketed Process Vessel 1000ml Vacuum Jacketed Process Vesse RR132002 2000ml Single Jacketed Process Vessel RR132052 2000ml Vacuum Jacketed Process Vesse RR133002 3000ml Single Jacketed Process Vessel RR133052 3000ml Vacuum Jacketed Process Vessel 5000ml Single Jacketed Process Vessel 5000ml Vacuum Jacketed Process Vessel PTFE Stirrer Shafts (for Process Vessels) Anchor PTFE Stirrer Shaft - 100ml Process Vessels Retreat Curve PTFE Stirrer Shaft - 100ml Process Vessels RR150107 Turbine PTFE Stirrer Shaft - 100ml Process Vessels Anchor PTFE Stirrer Shaft - 250ml Process Vessels Retreat Curve PTFE Stirrer Shaft - 250ml Process Vessels RR150257 Turbine PTFE Stirrer Shaft - 250ml Process Vessels RR150502 Anchor PTFE Stirrer Shaft - 500ml Process Vessels Retreat Curve PTFE Stirrer Shaft - 500ml Process Vessels RR150507 RR150512 RR151002 Turbine PTFE Stirrer Shaft - 500ml Process Vessels Anchor PTFE Stirrer Shaft - 1000ml Process Vessels Retreat Curve PTFE Stirrer Shaft - 1000ml Process Vessels RR151012 Turbine PTFE Stirrer Shaft - 1000ml Process Vessels Anchor PTFE Stirrer Shaft - 2000ml Process Vessels RR152002 RR152007 Retreat Curve PTFE Stirrer Shaft - 2000ml Process Vessels Turbine PTFE Stirrer Shaft - 2000ml Process Vessels RR152012 Anchor PTFE Stirrer Shaft - 3000ml Process Vessels Retreat Curve PTFE Stirrer Shaft - 3000ml Process Vessels Turbine PTFE Stirrer Shaft - 3000ml Process Vessels RR153007 RR153012 RR155002 Anchor PTFE Stirrer Shaft - 5000ml Process Vessels RR155007 Retreat Curve PTFE Stirrer Shaft - 5000ml Process Vessels Turbine PTFE Stirrer Shaft - 5000ml Process Vessels

AVA Lab Control Software - take control of your chemistry

Control multi-device jacketed reaction system

Multi-device reactions

- · Control up to 4 reaction systems with up to 16 devices on one screen
- Create complex experiments with any number of steps in series or in parallel
- Pre-program multi-step recipe, with the flexibility to make and track on-the-fly adjustments
- Interlink devices and set feedback/control loops, end-point conditions and safety limits

No more manual data logging or manual errors

- · Automatically record reaction parameters and log what you do, as you do it
- Repeat experiments accurately for reproducible and consistent results
- · Create reports in a few clicks or export data as a CSV file for further analysis
- Share results between users to improve research and collaboration

Safe unattended chemistry

- · Automate cooling during exothermic events
- Link devices such as balances and pumps for controlled, unattended reagent additions
- · Define safety overrides and cut-off conditions
- Configure audio and visual alarms

Data Hub

Using the Data Hub, AVA can control and log multiple 3rd party devices with RS232 serial interface.



Data Hub Specifications

- 2 x Pt100 temperature sensor ports with LEMO connection
- 1 x Ethernet port

AVA Software Kit includes:

- AVA Software
- Laptop controller
- 4 Port data hub

4 x RS232 ports

- 2 x Ethernet cables

- 1 year AVA Care Support

Download

Free AVA

Demo

www.radleys.com

Control a variety of devices

with multiple devices

Control and log reaction systems

Level 4

AVA software includes a library of pre-configured driver files allowing easy integration with a wide range of 3rd party devices.

Radleys also provide a 'New Driver Configuration Service' if required.



























Download and try the software for FREE

· Free software updates during support period, keeping software current

Learn how AVA software works - try before you buy

· Simulate control of devices

AVA Care Support

· Free support for 1st year

· Free priority email and telephone support · Free set up and application support

- · Set up apparatus and control experimental Schedules
- Share setups and schedules with other AVA users
- · Analyse results and create reports for real or simulated data
- Find out more about AVA Level 1-3 on our website

Accelerating chemistry

Radleys provide innovative chemistry equipment for safer, cleaner, greener and more productive chemical research.

Visit www.radleys.com to see our full range of chemistry productivity tools.



Shire Hill, Saffron Walden, Essex, CB11 3AZ, United Kingdom, t: +44 1799 513320 e: sales@radleys.com



