

## 24 Position Parallel Work-Up

Designed to make your parallel chemistry work-up and purification quick and easy.



## GreenHouse Work-Up Station™

Rapid, sequential, parallel purification in a 24 well, microtitre plate footprint using standard 3 ml or 6 ml filtration, phase separation and SPE columns

### Features

- Innovative, stacked column design allows sequential purification e.g. phase separation and SPE or filtration and SPE.
- Removable, upper support plate allows unit to be used for single or sequential purification.
- Clear gas cover allows pressure to be applied to increase flow rate of difficult samples. Accommodates all standard 3 ml and 6 ml tabbed or tabless SPE or filtration columns.
- Accepts all standard 24 well microtitre plates, as well as GreenHouse reaction blocks for fraction collection.
- Full range of filtration, phase separation and SPE columns including silica, SCX, aminopropyl and SAX acetate.

